Benchmark E



Prefilled Tubes and Bulk Beads: The Right Bead for Every Sample

- Wide range of kits available
 Variety of bead sizes and types to
 meet your sample needs
- TriplePure™ Quality
 Nuclease, Protease and Nucleic
 Acid free
- Prefilled 2.0ml tubes or bulk beads
 Convenience of prefilled tubes
 Bulk beads for filling your own

Benchmark Scientific's homogenizer bead kits come in a variety of materials and sizes for processing a wide range of samples. Prefilled for convenience, the unskirted 2.0ml tube kits will fit in all popular bead homogenizers including Benchmark Scientific's BeadBug™ and BeadBalaster™ 24. Tested to ensure they are RNase, DNase, protease, and nucleic acid free, the TriplePure Zirconium Bead Kits are supplied in five bead sizes to homogenize everything from bacteria to large tissues. Acid washed silica and stainless steel bead kits as well as specialty bead mixes and bulk beads are also available.

| Bead material | Bead Diameter | Prefilled 2.0ml tubes, pk/50 | Bulk cat. no./bottle | | | | | |
|--|---------------------------------------|------------------------------|----------------------|--|--|--|--|--|
| TriplePure Zirconium | 0.1mm | D1032-01 | D1132-01TP/250g | | | | | |
| TriplePure Zirconium | 0.5mm | D1032-05 | D1132-05TP/250g | | | | | |
| TriplePure Zirconium | 1.0mm | D1032-10 | D1132-10TP/250g | | | | | |
| TriplePure Zirconium | 1.5mm | D1032-15 | D1132-15TP/250g | | | | | |
| TriplePure Zirconium | 3.0mm | D1032-30 | D1132-30TP/250g | | | | | |
| TriplePure Zirconium | 0.1, 0.5, 1.0, 1.5, 3.0mm | D1032-SK* | NA | | | | | |
| Silica | 0.1mm | D1031-01 | D1131-01/200g | | | | | |
| Silica | 0.5mm | D1031-05 | D1131-05/200g | | | | | |
| Silica | 1.0mm | D1031-10 | D1131-10/200g | | | | | |
| Stainless Steel | 2.8mm | D1033-28 | D1133-28/1,000 balls | | | | | |
| Garnet, Zirconium | Garnet shards/6mm zirconium satellite | D1033-30G | NA | | | | | |
| Silica, Zirconium | Mixed size | D1034-MX | NA | | | | | |
| Empty tubes, with scre | of 50 | D1031-T20 | | | | | | |
| TriplePure Zirconium 0.1, 0.5, 1.0, 1.5, 3.0mm D1032-SK* NA Silica 0.1mm D1031-01 D1131-01/200g Silica 0.5mm D1031-05 D1131-05/200g Silica 1.0mm D1031-10 D1131-10/200g Stainless Steel 2.8mm D1033-28 D1133-28/1,000 Garnet, Zirconium Garnet shards/6mm zirconium satellite D1033-30G NA Silica, Zirconium Mixed size D1034-MX NA Empty tubes, with screw cap and sealing ring, nonskirted, pack of 50 Empty tubes, with screw cap and sealing ring, nonskirted, pack of 1,000 Empty tubes, with screw cap and sealing ring, nonskirted, pack of 1,000 BeadBug Bead Homogenizer, three tube capacity D1030 | | | | | | | | |
| BeadBug Bead Homogo | | D1030 | | | | | | |
| Silica, Zirconium Mixed size D1034-MX NA Empty tubes, with screw cap and sealing ring, nonskirted, pack of 50 D1031-T20 Empty tubes, with screw cap and sealing ring, nonskirted, pack of 1,000 D1031-T21 | | | | | | | | |

| | Bead material | Silica | | | Zirconium | | | | Steel | Garnet with Zr Satellite | Silica and Zirconium | |
|----------------------|---------------------|--------|-------|-------|-----------|-------|-------|-------|-------|-----------------------------|-------------------------|-------|
| | Bead diameter | 0.1mm | 0.5mm | 1.0mm | 0.1mm | 0.5mm | 1.0mm | 1.5mm | 3.0mm | 2.8mm | 0.5/6mm | Mixed |
| | Adipose | | | | | | | Х | | | | |
| | Artery | | | | | | | Х | X | | | |
| | Bone | | | | | | | | | X | | |
| | Brain | | | | | | | X | Х | | | |
| | Insect - fly | | | | | | | | X | Х | | |
| | Insect - Tick | | | | | | | | Х | Х | | |
| | Insect - Mosquito | | | | | | X | X | X | | | |
| Animal | Gonad | | | | | | | Х | Х | | | |
| | Hair | | | | | | | | X | X | Х | |
| | Heart | | | | | | | Х | X | | | |
| | Kidney | | | | | | | Х | | | | |
| | Liver | | | | | | | | | | | |
| | Lung | | | | | | | X | X | | | |
| | Muscle | | | | | | | | Х | Х | | |
| | Pancreas | | | | | | | X | X | | | |
| | Skin | | | | | | | | X | | Х | |
| | Spleen | | | | | | | X | X | | | |
| | Teeth | | | | | | | | X | Х | | |
| Microbial - | Algae | | Х | | | Х | | | | | | |
| | Bacteria | X | Х | | Х | Х | | | | | | |
| | Fungi | | X | Х | | Х | Х | | | | | |
| | Mold | | X | | | Х | | | | | | |
| | Spores | | X | Х | | Х | Х | | | | | |
| | Yeast | | X | Х | | Х | Х | | | | | |
| Plant and Soil | Leaves | | | | | | | | Х | Х | | |
| | Nuts | | | | | | | | | Х | | |
| | Plant, soft tissue | | | | | | | | X | Х | | |
| | Plant, tough tissue | | | | | | | | X | Х | Х | |
| | Seeds | | | | | | | | X | Х | | |
| | Stems | | | | | | | | X | X | | |
| | Soil | | | | | | | Х | | | | Х |
| | Coral | | | | | | | | | X | | |
| Fecal | Stool | | | | | | | X | | | | Х |
| | Scat | | | | | | | | Х | Х | | Х |

Above is a partial list of sample types and suggested beads. For additional information and technical support, contact Benchmark Scientific at info@BenchmarkScientific.com.



PO Box 709 Edison, NJ 08818 TEL: 908 769-5555 Fax: 908 222-1864 email: info@BenchmarkScietific.com Web: www.BenchmarkScientific.com